

TECHNICAL DATA SHEET

Part#: 5-S-SE



Tweezers - 1 Star Straight Tapered Ultra Fine Point - SS

The 5-S-SE straight tapered ultra-fine pointed tweezers is primarily used for microscopy. The style 5 has the finest points available. The overall length is 4.25". It is made of stainless steel.

• Precision Level	★
• Material	SS (AISI 440)
• Surface Resistivity	Conductive
• Shape	Straight Tapered
• Group	Conventional Style Tweezers
• Weight (OZ)	1
• Tip Size (CTM)	N/A
• Tip Angle	Straight
• Tip Length	
• Body Width	.375" (9.38mm)
• Body Height One Leg	.06" (1.5mm)
• OAL	4.25" (106.25mm)
• Country of Origin	Pakistan
• UPC	662847002700
• Schedule B	8203.20.20.00
• Harmonization Code	8203.20.20
• ECCN	EAR99
• UNSPSC Code	27112105



C.T.M.
Critical Tip Measurements
(Tip Width and Thickness of one point)



EXCELTA[®]
CORPORATION

A) N/A
B) N/A
C) 0.375"
D) 4.25"
E) 0.060"