

TECHNICAL DATA SHEET

Part#: 5-SA



Tweezers - 3 Star Straight Tapered Ultra Fine Point - Anti-Mag. SS

The 5-SA straight tapered ultra-fine pointed tweezers is primarily used for microscopy. The style 5 has the finest points available. The overall length is 4.25". It is made of anti-magnetic Neverust[®] stainless steel.

- Precision Level ★★★
- Material Anti-Mag. SS (316L)
- Surface Resistivity Conductive
- Shape Straight Tapered
- Group Conventional Style Tweezers
- Weight (OZ) 1
- Tip Size (CTM) .004" x .002" (.1mm x .05mm)
- Tip Angle Straight
- Tip Length
- Body Width .375" (9.38mm)
- Body Height One Leg .06" (1.5mm)
- OAL 4.25" (106.25mm)
- Country of Origin Switzerland
- UPC 662847002717
- Schedule B 8203.20.20.00
- Harmonization Code 8203.20.20
- ECCN EAR99
- UNSPSC Code 27112105



C.T.M.
Critical Tip Measurements
(Tip Width and Thickness of one point)



EXCELTA[®]
CORPORATION

- 5-SA ★★★
- A) 0.0039"
 - B) 0.0024"
 - C) 0.375"
 - D) 4.25"
 - E) 0.060"