

## **TECHNICAL DATA SHEET**

Part#: AF-1-SE

Forceps - Laboratory - Straight .03" x .05" Serrated Tips - SS

This Adson tissue 4.75" forceps is used for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The serrated tips provide more 'sure' grip. The satin finish reduces glare.

Precision Level

• Material SS (AISI 440)

Surface ResistivityShapeConductiveStraight Tapered

• Group Tweezers (Forceps)

• Weight (OZ)

• Tip Size (CTM) .03" x .05" (.75mm x 1.25mm)

• Tip Angle Straight

Tip Length

Body Width 0.43"(1.08mm)Body Height One Leg 0.04"(1mm)

• OAL 4.75" (118.75mm)

• Country of Origin Pakistan

• UPC 662847011924 • Schedule B 8203.20.20.00

Harmonization Code 8203.20.20ECCN EAR99

• UNSPSC Code 27112105



