

TECHNICAL DATA SHEET

Part#: AF-3-SE

Forceps - Laboratory - Straight .083" x .074" Serrated Tips - SS

This Adson tissue 4.75" large forceps (1 X 2 teeth) is used for holding and manipulating delicate tissues. These forceps have a wide, flat thumb grasp area that is commonly serrated. The satin finish reduces glare.

Precision Level

• Material SS (AISI 440)

Surface Resistivity
 Shape
 Group
 Conductive
 Straight Tapered
 Tweezers (Forceps)

Weight (OZ)

• Tip Size (CTM) .083" x .074" (2.08mm x 1.85mm

27112105

• Tip Angle Straight

Tip Length

• UNSPSC Code

• Body Width 0.43"(1.08mm)
• Body Height One Leg 0.04"(1mm)

• OAL 4.75" (118.75mm)

Country of Origin
 UPC
 Schedule B
 Harmonization Code
 ECCN
 Pakistan
 662847011931
 8203.20.20.00
 8203.20.20
 EAR99



