

TECHNICAL DATA SHEET

Part#: F-7-SE

Forceps - Laboratory - Straight .02" x .03" Serrated Tips - SS

The F-7-SE straight 4.25" dressing forceps are used for holding dressing materials during procedures. It has serrated tips for better gripping.

• Precision Level	★
• Material	SS (AISI 440)
• Surface Resistivity	Conductive
• Shape	Straight
• Group	Tweezers (Forceps)
• Weight (OZ)	0.25
• Tip Size (CTM)	.02" x .030" (.5mm x .75mm)
• Tip Angle	Straight
• Tip Length	
• Body Width	0.385"(9.63mm)
• Body Height One Leg	0.063"(1.58mm)
• OAL	4.25" (106.25mm)
• Country of Origin	Pakistan
• UPC	662847012655
• Schedule B	8203.20.20.00
• Harmonization Code	8203.20.20
• ECCN	EAR99
• UNSPSC Code	27112105

